TEST NAME/ PSYCHE CODE	COLLECTION GUIDELINES	REFERENCE VALUES	ROOM TEMP.	REFRIGERATED	FROZEN	Clinical Utility
ALBUMIN, SERUM	2.0 ml SERUM. 8-12HR FAST	3.5-5.5 g/dL		48 hours		To screen for a liver disorder or kidney disease or to evaluate nutritional status along with or
CODE: ALB	PREFERRED.					instead of a prealbumin test.
CPT: 82040	SCHEDULE: DAILY					
ALCOHOL, URINE	30 ml RANDOM URINE	NEGATIVE		6-24 hours		To test for a certain level of Alcohol in the system
CODE: UALC	30 IIII KANDOWI OKINE	NEGATIVE				
CPT: 82055	SCHEDULE: DAILY					
ALKALINE PHOSPHATASE	1 ml UNHEMOLYZED SERUM.	38-126 IU/L		48 hours		To screen for or monitor treatment of a liver or bone disorder.
CODE: ALKP	THII GIVITEWOLIZED SERGIVI.	36-120 10/1				
CPT:84075	SCHEDULE: DAILY					
ALT, SERUM (SGPT)	1 ml UNHEMOLYZED SERUM.	FEMALE: 4-48 IU/L		48 hours		An alanine aminotransferase (ALT) test measures the amount of this enzyme in theblood.
CODE: ALT	THI GIVILLINGETZED SERGIVI.	MALE: 8-36 IU/L				Increased in viral hepatitis and myocardial infarction,
CPT: 84460	SCHEDULE: DAILY					
AMMONIA (NH3)	TESTING PERFORMED WITH STRICT		Green Top-		Lavendar-7 days	Measures the amount of ammonia in the blood, formed by decoposition of urea.
	GUIDELINES: MUST BE TO LAB		must be put on			
	WITHIN 30 MIN OF COLLECTION.		ice 30 mins after			
	GREEN TOP (HEPARANIZED) WHOLE		collection			_
	BLOOD, KEPT IN A SLURRY OF					
	WATER AND ICE.					
	IF GREATER THAN 30 MIN: COLLECT	11-35 umol/L (IN-HOUSE TEST)				
	LAVENDER EDTA WHOLE BLOOD.	<48 umol/L (SENDOUT TEST)				
	SPIN, POUR OFF EDTA PLAMSA AND					
CODE: NH3	FREEZE WITHIN 20 MINUTES OF COLLECTION (TO BE SENT TO					
	REFERENCE LAB.)					
	PATIENT SHOULD NOT UNDERGO					
	VIGOROUS MUSCULAR ACTIVITY					
	PRIOR TO TESTING.					
CPT: 82140	SCHEDULE: DAILY					
CF1: 82140	SCHEDULE: DAILY		<u> </u>		<u> </u>	To measure the amount of AMPHETAMINE in the blood.
AMPHETAMINE, URINE QUALITATIVE	30 ml RANDOM URINE	INTERPRETIVE REPORT		2-4 days		To measure the amount of Alvient TAIVIIIVE III the blood.
CODE: AMP				,.		
CPT: 82145	SCHEDULE: DAILY					
AMYLASE, SERUM	1 0 ml UNHEMOLYZED SERUM	ml unhemolyzed serum 30-70 U/L		48 hours		An amylase test measures the amount of this enzyme which is necessary for starch digestion.
CODE: AMYL	1.5 III GIVILLINGETZED SEKGIVI					
CPT: 82150	SCHEDULE: DAILY					
ASPIRIN RESISTANCE	2 ml 3.2% NA CITRATE (PEDI BLUE	GREATER THAN 550= ASPIRIN				Measures the effectiveness of aspirin therapy. Some people are 'non-responsive' to aspirin therapy and require treatment with anotehr ancticoagulant.
	TOP) SPECIMEN MUST BE RUN WITHIN 4 HOURS OF COLLECTION.	RESISTANT	STANT			
CODE: BPA		LESS THAN 550= ASPIRIN NON-				
CPT: 85576	SCHEDULE: DAILY	RFSISTANT				
ASPARTATE AMINOTRANSFERASE	1 ml UNHEMOLYZED SERUM	7-40 IU/L		48 hours		An aspartate aminotransferase (AST) test measures the amount of this enzyme in the blood. AST
AST, SERUM (SGOT)						is normally found in red blood cells, liver, heart, muscle tissue, pancreas, and kidneys. Elevated in
CODE: AST						viral hepatitis and mycardial infarction.
CPT: 84450	SCHEDULE: DAILY					