TEST NAME/ PSYCHE CODE	COLLECTION GUIDELINES	REFERENCE VALUES	ROOM TEMP.	REFRIGERATED	FROZEN	Clinical Utility
INSULIN	1 ml SERUM	1.9-23.0 IU/ml	8 hours	24 hours	UNDETERMINED	Insulin may be used, along with glucose and C-peptide levels, to help diagnose insulinomas and to help diagnose the cause of documented acute or chronic hypoglycemia.
CODE: INSU						
CPT: 83525	SCHEDULE: DAILY					
IRON TOTAL	1 mal CEDUM	25-170 ug/dL		48 hours		This test helps determine your blood iron level. Iron is increased in chronic iron overload diseases,
CODE: IRON	1 ml SERUM					particularly hereditary hemochromatosis.
CPT: 83540	SCHEDULE: DAILY					
IRON PANEL (IRON TOTAL AND BINDING CAPACITY)	2 ml SERUM	IRON=25-170 ug/dL TIBC=250-350 ug/dL SATURATION= 12-57%		48 hours		Screens for iron deficiency
CODE: IRNPNL		Transferrin				
CPT: 7573	SCHEDULE: DAILY					