

TEST NAME/ PSYCHE CODE	COLLECTION GUIDELINES	REFERENCE VALUES	ROOM TEMP.	REFRIGERATED	FROZEN	Clinical Utility
OCCULT BLOOD (HEMOCCULT, FECAL OCCULT BLOOD)  CPT: 82272	STOOL SPECIMEN OR STOOL CARD	NEGATIVE	24 hours	-----	-----	The main use for the fecal occult blood test (FOBT) or the fecal immunochemical test (FIT) is as a screen for early colon cancer.
	Do not refridgerate or freeze SCHEDULE: DAILY					
OPIATES, URINE QUALITATIVE  CODE: OPIQ CPT: 80101	30 ml URINE	INTERPRETIVE REPORT	-----	2-4 days	-----	Testing detects this substance not normally found in the body
	SCHEDULE: DAILY					