

TEST NAME/ PSYCHE CODE	COLLECTION GUIDELINES	REFERENCE VALUES	ROOM TEMP.	REFRIGERATED	FROZEN	Clinical Utility
FECAL LEUKOCYTES CODE: LEUF CPT: 85540	1 ml FRESH STOOL  SCHEDULE: DAILY	NONE SEEN	-----	24 hours	-----	Used to detect WBC's in stools indicating infection and colitis
FERRITIN CODE: FERR CPT: 82728	2 ml UNHEMOLYZED SERUM.  SCHEDULE: DAILY	MALE: 23.9-336.2 ng/mL FEMALE: 11.0-306.8 ng/mL	8 hours	48 hours	UNDETERMINED	The ferritin test is ordered to assess a person's iron stores in the body
FOLATE (FOLIC ACID), SERUM CODE: FOLT  CPT: 82746	1 ml UNHEMOLYZED SERUM. SUBMIT REFRIDGERATED OR FROZEN.  SCHEDULED: DAILY	3.0-20.0 ng/mL	Refig. ASAP	8 hours	6 months	Folate is primarily ordered to detect deficiencies and to help diagnose the cause of certain anemias
FREE T3 CODE: FT3 CPT: 84481	1 ml SERUM  SCHEDULE: DAILY	2.5-3.9 pg/mL	8 hours	24 hours	UNDETERMINED	A T3 test is used to assess thyroid function. It is ordered primarily to help diagnose hyperthyroidism and may be ordered to help monitor the status of a person with a known thyroid disorder.
FREE TESTOSTERONE  CODE: TFTS  CPT: 37073	1.0 ML SERUM  SCHEDULE: DAILY	MALE: 20-50YR 24.3-11.02% FEMALE: 20-46YR 065-10.93% POST MENOPAUSE: 47-97YR .23-6.80%	8 hours	7 days	7 days	Results indicate hormone replacement adequacy or necessity. Decreased levels may occur in liver disease or severe systemic illness.
FSH, SERUM (FOLLICLE STIMULATING HORMONE) CODE: FSH CPT: 83001	1 ml SERUM. REFRIGERATE  SCHEDULE: DAILY	MALE: 1.0-12.0 mIU/mL FEMALE: 3.0-20.0 mIU/mL	8 hours	48 hours	3 months	To evaluate the pituitary function, especially as it relates to fertility issues, gonadal failure, maturation concerns, or pituitary tumors
FUNGAL STAINS CODE: FUNST CPT: 87206	  SCHEDULE: DAILY	NO FUNGAL FORMS SEEN	-----	48 HOURS	-----	Rapid detect of fungal forms in specimens.